

BIL	Quotation Reference	Description	Advertisement Date	Closing Date (Not Later Than 2.00PM)	Quotation Fee	Requesting Department	Focal Person
1	BME/A2K/(930)2022/HTD(A4)	<p style="text-align: center;"><b>TO REPAIR ONE (1) UNIT OF ULTRASOUND MACHINE, BRAND: SONOSITE, MODEL: S-NERVE AT OPERATING THEATRE</b></p> <p style="text-align: center;"><u>PLACE OF SUBMISSION:</u> QUOTATION BOX (GROUND FLOOR) MINISTRY OF HEALTH COMMONWEALTH DRIVE BANDAR SERI BEGAWAN, BB 3910 NEGARA BRUNEI DARUSSALAM</p>	25/10/2022	12/11/2022	\$5.00	JABATAN TEKNOLOGI PENJAGAAN KESIHATAN, KEMENTERIAN KESIHATAN.	<p style="text-align: center;">MUHD AMIRUL FAZLEEN BIN HAJI KHALIDIN</p> <p style="text-align: center;">JURUTERA PERUBATAN HAYAT</p> <p style="text-align: center;">TEL: 2242424 ext. 8693 (RIPASH)</p>

# LAMPIRAN - A1

RUJ: BME/A2K/(930)/2022/HTD (B4)

11/10/2022

TO REPAIR ONE (1) UNIT OF ULTRASOUND MACHINE, BRAND : SONOSITE, MODEL : S-NERVE AT OPERATING THEATRE

## SCOPE OF WORK

### EQUIPMENT DETAILS:

EQUIPMENT: ULTRASOUND

BRAND: SONOSITE

MODEL: S-NERVE

BME NO : 7030, SN : 03G2FI

LOCATION : OPERATING THEATRE

### SPAREPART DETAILS:

I) BATTERY PACK [P/N:35770P] - 1 UNIT

II) LOW NOISE POWER SUPPLY [P/N:32890P]- 1 UNIT

### SCOPE OF WORKS

- 1) DISSEMBLE THE UNIT FOR FURTHER INSPECTION
- 2) CLEAN THE INTERIOR MECHANICAL PARTS OF THE UNIT
- 3) REMOVE THE FAULTY AND BROKEN POWER SUPPLY OF THE UNIT
- 4) REPLACE THE POWER CABLE
- 5) REASSEMBLE THE UNIT
- 6) TEST THE UNIT WITH THE NEWLY INSTALLED POWER SUPPLY

## LAMPIRAN - A2

RUJ: BME/A2K/(930)/2022/HTD (B4)

TO REPAIR ONE (1) UNIT OF ULTRASOUND MACHINE, BRAND : SONOSITE, MODEL : S-NERVE AT OPERATING THEATRE

NO	SPAREPART	PART NUMBER	QUANTITY
1	BATTERY PACK	35770P	1
2	LOW NOISE POWER SUPPLY	32890P	1