

BIL	Quotation Reference	Description	Advertisement Date	Closing Date (Not Later Than 02.00PM)	Quotation Fee	Requesting Department	Focal Person
3	BME/A2K / (1358)/2025/HTD (A20)	<p>ONE OFF PREVENTIVE MAINTENANCE FOR ONE (1) UNIT OF PROTENI / NITROGEN ANALYZER, BRAND: FOSS, MODEL: KJELTEC 8400, S/N: 9188846, BME2905 AT FOOD CHEMISTRY LABORATORY, DEPARTMENT OF SCIENTIFIC SERVICES</p> <p>PLACE OF SUBMISSION: QUOTATION BOX (GROUND FLOOR) MINISTRY OF HEALTH COMMONWEALTH DRIVE BANDAR SERI BEGAWAN, BB 3910 NEGARA BRUNEI DARUSSALAM</p>	05/02/2025	22/02/2025	\$5.00	<p>JABATAN TEKNOLOGI PENJAGAAN KESIHATAN, KEMENTERIAN KESIHATAN.</p>	<p>HAJI ZAENUL 'ARIFIEN BIN DATO PADUKA HAJI MOHD ROSLI</p> <p>JURUTERA PERUBATAN HAYAT</p> <p>TEL: 2382424 ext. 6122 / 6232</p>

LAMPIRAN - A1

RUJ: BME/A2K/(1358)/2025/HTD (A20)

22/1/2025

ONE OFF PREVENTIVE MAINTENANCE FOR ONE (1) UNIT OF PROTEIN/ NITROGEN ANALYZER, BRAND: FOSS, MODÈL: KJELTEC 8400, S/N: 91888846, BME29405 AT FOOD CHEMISTRY LABORATORY, DEPARTMENT OF SCIENTIFIC SERVICES

SCOPE OF WORK

One-off preventive maintenance services for one (1) unit of Protein/ Nitrogen Analyzer System, Brand: Foss, Model: Kjelttec 8400, S/N: 91888846, BME29405, at Food Chemistry Laboratory, Department of Scientific Services, Berakas Branch, inclusive of but not limited to:

The system includes the following instruments:

1. Kjelttec 8400 System, S/N: 91888846
2. Kjelttec 8420 Autosampler, S/N: 91889497
3. Digestor 2520 System, S/N: 91888097
4. Scrubber 2501 System, S/N: 91889504

- Scope of work to follow manufacturer's preventive maintenance manual/recommendation specific for the respective instruments listed above, which include the following where necessary:

- (i) inspection
- (ii) cleaning
- (iii) installation of PM kits, consumables and spare parts listed in spareparts list (Lampiran A2)
- (iv) alignment
- (v) calibration
- (vi) any other related preventive maintenance works required.

- Functional check post-maintenance.

- Apply preventive maintenance sticker/tag after successful service.

- Disposal of any used spare parts/pm kits, unless otherwise instructed.

- Work to be carried out by engineer(s) authorised in providing maintenance to the equipment.

- Vendor(s) to submit maintenance checklist specific for the equipment along with other necessary documents to join this tender.

- Upon approval, the earliest the preventive maintenance work can be carried out on the equipment is one month before the next PM due date, and not anytime earlier.

LAMPIRAN - A2

RUJ:

BME/A2K/(1358)/2025/HTD (A20)

ONE OFF PREVENTIVE MAINTENANCE FOR ONE (1) UNIT OF PROTEIN/ NITROGEN ANALYZER, BRAND: FOSS, MODEL: KJELTEC 8400, S/N: 91888846, BME29405 AT FOOD CHEMISTRY LABORATORY, DEPARTMENT OF SCIENTIFIC SERVICES

	SPAREPART	PART NUMBER	QUANTITY
1	FOSS PM Kit, KT8400		1 set
2	FOSS PM Kit, KT8420		1 set
3	Grinding head for Digestor		1 set
4	FOSS PM Kit, Scrubber		1 set
5	Temperature Cut-Out DS20		1 set