Rujukan Kami: (86) MOH/HQ/P/IKLAN-SH/2025

LAMPIRAN 1

BIL	Quotation Reference	Description	Advertisement Date	Closing Date (Not Later Than 09.00AM)	Quotation Fee	Requesting Department
1	BME/A2K/(1443)/2025/HTD (A17)	TO SERVICE ONE OFF PREVENTIVE MAINTENANCE (PPM) FOR ONE (1) UNIT OF AUTOMATED HIV ANALYZER, BRAND: INNIGENETICS, MODEL: AUTO-LIA 48 INCLUDING ACER XC1785 COMPUTER (WITH INSTALLATION OF THE LATEST SOFTWARE LIRAS V3.00) AND AN ADVANCE A4 FLATBED SCANNER. (WO: 2024- 4313)	10/06/2025	05/07/2025	\$5.00	JABATAN TEKNOLOGI PENJAGAAN KESIHATAN, KEMENTERIAN KESIHATAN

28/5/2025

RUJ: BME/A2K/(1443)/2025/HTD (A17)

TO SERVICE ONE OFF PREVENTIVE MAINTENANCE (PPM) FOR ONE (1) UNIT OF AUTOMATED HIV ANALYZER, BRAND: INNOGENETICS, MODEL: AUTO-LIA 48 INCLUDING ACER XC1785 COMPUTER (WITH INSTALLATION OF THE LATEST SOFTWARE LIRAS V3.00) AND AN ADVANCED A4 FLATBED SCANNER. (WO: 2024-4313)

SCOPE OF WORK

EQUIPMENT: AUTOMATED HIV ANALYZER

BRAND: INNOGENETICS MODEL: AUTO-LIA 48

S/N: 501000002 BME NO: BME5363

WORK ORDER: 2024-4313 LOCATION: SUMBILING LAB.

PARTS TO BE REPLACED:

DISPENSING TUBES CALIBRATION TUBING SCOPE OF WORK:

- Verify System using floppy diskette of the system
- Check & clean optional printer.
- Check & clean Display Monitor Module
- Clean & grase guiding bars
- Check & clean pump head
- Clean verifier module temperature
- Clean & perform pump test
- Change calibration tubing
- Clean level sensor pin intubing
- Check & clean front cooling fan
- Check AC/DC power supply voltage
- Replace all dispensing tubes
- Perform Tub volume adjustment
- Perform battery measurement
- Perform pump calibration
- Perform pump channel adjustment
- Perform pump volume recalibration
- General cleaning on the instrument

LAMPIRAN - A2

RUJ: BME/A2K/(1443)/2025/HTD (A17)

TO SERVICE ONE OFF PREVENTIVE MAINTENANCE (PPM) FOR ONE (1) UNIT OF AUTOMATED HIV ANALYZER, BRAND: INNOGENETICS, MODEL: AUTO-LIA 48 INCLUDING ACER XC1785 COMPUTER (WITH INSTALLATION OF THE LATEST SOFTWARE LIRAS V3.00) AND AN ADVANCED A4 FLATBED SCANNER. (WO: 2024-4313)

	SPAREPART	PART NUMBER	QUANTITY
1	DISPENSING TUBES		1
2	CALIBRATION TUBING		1