

QTN REF: (28) PPN/QTN/2025 (POCD)

**SUPPLY AND DELIVERY OF RESTORATIVE MATERIALS AND INSTRUMENTS FOR COMPOSITE RESTORATION PROCEDURES FOR DEPARTMENT OF DENTAL SERVICES,
MINISTRY OF HEALTH**

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
1.	<p>High Speed White Stone Burs (Flame, FL2)</p> <p>Specifications:</p> <ul style="list-style-type: none"> • Flame shape • Made of micrograined alumina oxide grits • For finishing of enamel, glass ionomer cements, composites, porcelain and prepared tooth abutment • Friction grip • ISO 025 (autoclavable) • Sample must be submitted upon participation for accurate evaluation 	70 Boxes (Boxes of 12s)				
2.	<p>High Speed White Stone Burs (Cone, CN1)</p> <p>Specifications:</p> <ul style="list-style-type: none"> • Cone shape • Made of micrograined alumina oxide grits • For finishing of enamel, glass ionomer cements, composites, porcelain and prepared tooth abutment • Friction grip • ISO 024 (autoclavable) • Sample must be submitted upon participation for accurate evaluation 	70 Boxes (Boxes of 12s)				

QTN REF: (28) PPN/QTN/2025 (POCD)

**SUPPLY AND DELIVERY OF RESTORATIVE MATERIALS AND INSTRUMENTS FOR COMPOSITE RESTORATION PROCEDURES FOR DEPARTMENT OF DENTAL SERVICES,
MINISTRY OF HEALTH**

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
3.	<p>High Speed White Stone Burs (Round, RD1)</p>  <p>Specifications:</p> <ul style="list-style-type: none"> • Round shape • Made of micrograined alumina oxide grits • For finishing of enamel, glass ionomer cements, composites, porcelain and prepared tooth abutment • Friction grip • ISO 030 (autoclavable) • Sample must be submitted upon participation for accurate evaluation 	70 Boxes (Boxes of 12s)				
4.	<p>Nanohybrid Composite Shade B3</p> <p>Specifications:</p> <ul style="list-style-type: none"> • Use in all classes of direct restorations (anterior or posterior), core build ups, splinting and indirect restorations (onlays, inlays and veneers) • Easy to polish for predictable esthetics (retains its natural gloss and produce natural-looking results) • Great handling without slumping and sticking to instruments • Holds its shape prior to cure • High strength with long lasting wear resistance • Delivery system is in capsular form (refills) • Long shelf life (expiry date should be at least 18 months from delivery date) • Sample must be submitted upon participation for accurate evaluation 	50 Packs (Of 20 Compoules)				

SUPPLY AND DELIVERY OF RESTORATIVE MATERIALS AND INSTRUMENTS FOR COMPOSITE RESTORATION PROCEDURES FOR DEPARTMENT OF DENTAL SERVICES,
MINISTRY OF HEALTH

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
5	<p>Nanohybrid Composite Shade C3</p> <p>Specifications:</p> <ul style="list-style-type: none"> • Use in all classes of direct restorations (anterior or posterior), core build ups, splinting and indirect restorations (onlays, inlays and veneers) • Easy to polish for predictable esthetics (retains its natural gloss and produce natural-looking results) • Great handling without slumping and sticking to instruments • Holds its shape prior to cure • High strength with long lasting wear resistance • Delivery system is in capsular form (refills) • Long shelf life (expiry date should be at least 18 months from delivery date) • Sample must be submitted upon participation for accurate evaluation 	<p>50 Packs (Of 20 Compoules)</p>				
6	<p>Clear Matrix Strip (6mm)</p> <ul style="list-style-type: none"> • Transparent strips used for composite restorations. • Manufactured from mylar polyester film • Provides greater strength than acetate • Resistant to solvents • Does not stretch, curl or break • Each strip with a thickness of 0.05mm • Sample must be submitted upon participation for accurate evaluation 	<p>60 Boxes</p>				

QTN REF: (28) PPN/QTN/2025 (POCD)

**SUPPLY AND DELIVERY OF RESTORATIVE MATERIALS AND INSTRUMENTS FOR COMPOSITE RESTORATION PROCEDURES FOR DEPARTMENT OF DENTAL SERVICES,
MINISTRY OF HEALTH**

NO	ITEM DESCRIPTION	QUANTITY	BRAND	MANUFACTURER AND ORIGIN	UNIT PRICE	TOTAL PRICE
7	<p>Clear Matrix Strip (8mm)</p> <p>Specifications:</p> <ul style="list-style-type: none"> • Transparent strips used for composite restorations. • Manufactured from mylar polyester film • Provides greater strength than acetate • Resistant to solvents • Does not stretch, curl or break • Each strip with a thickness of 0.05mm. • Sample must be submitted upon participation for accurate evaluation 	60 Boxes				
8	<p>Clear Matrix Strip (10 Mm)</p> <p>Specifications:</p> <ul style="list-style-type: none"> • Transparent strips used for composite restorations. • Manufactured from mylar polyester film • Provides greater strength than acetate • Resistant to solvents • Does not stretch, curl or break • Each strip with a thickness of 0.05mm. • Sample must be submitted upon participation for accurate evaluation 	60 Boxes				

QTN REF: (28) PPN/QTN/2025 (POCD)

SUPPLY AND DELIVERY OF RESTORATIVE MATERIALS AND INSTRUMENTS FOR COMPOSITE RESTORATION PROCEDURES FOR DEPARTMENT OF DENTAL SERVICES,
MINISTRY OF HEALTH

DELIVERY OF PERIOD: Not more than 60 days upon confirmation	(Yes/No) (If No, please specify)
<p>PRICE VALIDITY:</p> <p><i>The quotation shall remain valid for 6 MONTHS from the final date for the submission of the quotation and no supplier may withdraw his/her quotation with that period. The Government reserves the right to extend this period if deemed necessary provided that such extension to the quotation validity period shall have the written consent of the suppliers.</i></p>	

QTN REF: (28) PPN/QTN/2025 (POCD)

SUPPLY AND DELIVERY OF RESTORATIVE MATERIALS AND INSTRUMENTS FOR COMPOSITE RESTORATION PROCEDURES FOR DEPARTMENT OF DENTAL SERVICES,
MINISTRY OF HEALTH

TERMS AND CONDITIONS			
a.	Tenderer must be registered with the Ministry of Health	<i>Acknowledgement:</i> <i>Company Ref No:</i>	Company's Official Stamp
b.	Please fill in the QUOTATION FORM completely including the USER REQUIREMENT FORM (if available). Submission of incomplete form may cause DISQUALIFICATION OF QUOTATION	I hereby certify the above quote to be correct. Signature: Name: 	
c.	Each tenderer is allowed to quote ONE BRAND WITH ONE PRICE ONLY for each item, Submission of more than one brand and price will cause DISQUALIFICATION OF QUOTATION.	Designation: 	
d.	Delivery Period: NOT MORE THAN 60 DAYS	Date: 	
e.	Please do not use TIPPEX for amendment.	

**QTN REF: (28) PPN/QTN/2025 (POCD)
 SUPPLY AND DELIVERY OF RESTORATIVE MATERIALS AND INSTRUMENTS FOR COMPOSITE RESTORATION
 PROCEDURES FOR DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH**

NO	ITEM DESCRIPTION	USER REQUIREMENT	YES	NO	PLEASE SPECIFY THE SPECIFICATION / SCOPE OFFERED FAILURE TO FILL THE OFFER DETAIL SHALL BE IDENTIFIED AS NON-COMPLIANCE TO THE OFFER
1.	High Speed White Stone Burs (Flame, FL2)	Flame shape			
		Made of micrograined alumina oxide grits			
		For finishing of enamel, glass ionomer cements, composites, porcelain and prepared tooth abutment			
		Friction grip			
		ISO 025 (autoclavable)			
		Sample must be submitted upon participation for accurate evaluation			
		70 Boxes (Boxes of 12s)			
2.	High Speed White Stone Burs (Cone, CN1)	Cone shape			
		Made of micrograined alumina oxide grits			
		For finishing of enamel, glass ionomer cements, composites, porcelain and prepared tooth abutment			
		Friction grip			
		ISO 024 (autoclavable)			
		Sample must be submitted upon participation for accurate evaluation			
		70 Boxes (Boxes of 12s)			

QTN REF: (28) PPN/QTN/2025 (POCD)
SUPPLY AND DELIVERY OF RESTORATIVE MATERIALS AND INSTRUMENTS FOR COMPOSITE RESTORATION
PROCEDURES FOR DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH

NO	ITEM DESCRIPTION	USER REQUIREMENT	YES	NO	PLEASE SPECIFY THE SPECIFICATION / SCOPE OFFERED FAILURE TO FILL THE OFFER DETAIL SHALL BE IDENTIFIED AS NON-COMPLIANCE TO THE OFFER
3.	High Speed White Stone Burs (Round, RD1)	Round shape			
		Made of micrograined alumina oxide grits			
		For finishing of enamel, glass ionomer cements, composites, porcelain and prepared tooth abutment			
		Friction grip			
		ISO 030 (autoclavable)			
		Sample must be submitted upon participation for accurate evaluation			
		70 Boxes (Boxes of 12s)			
4.	Nanohybrid Composite Shade B3	Use in all classes of direct restorations (anterior or posterior), core build ups, splinting and indirect restorations (onlays, inlays and veneers)			
		Easy to polish for predictable esthetics (retains its natural gloss and produce natural-looking results)			
		Great handling without slumping and sticking to instruments			
		Holds its shape prior to cure			
		High strength with long lasting wear resistance			
		Delivery system is in capsular form (refills)			
		Long shelf life (expiry date should be at least 18 months from delivery date)			
		Sample must be submitted upon participation for accurate evaluation			
		Use in all classes of direct restorations (anterior or posterior), core build ups, splinting and indirect restorations (onlays, inlays and veneers)			
50 Packs (Of 20 Compoules)					

QTN REF: (28) PPN/QTN/2025 (POCD)
SUPPLY AND DELIVERY OF RESTORATIVE MATERIALS AND INSTRUMENTS FOR COMPOSITE RESTORATION
PROCEDURES FOR DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH

NO	ITEM DESCRIPTION	USER REQUIREMENT	YES	NO	PLEASE SPECIFY THE SPECIFICATION / SCOPE OFFERED FAILURE TO FILL THE OFFER DETAIL SHALL BE IDENTIFIED AS NON-COMPLIANCE TO THE OFFER
5.	Nanohybrid Composite Shade C3	Use in all classes of direct restorations (anterior or posterior), core build ups, splinting and indirect restorations (onlays, inlays and veneers)			
		Easy to polish for predictable esthetics (retains its natural gloss and produce natural-looking results)			
		Great handling without slumping and sticking to instruments			
		Holds its shape prior to cure			
		High strength with long lasting wear resistance			
		Delivery system is in capsular form (refills)			
		Long shelf life (expiry date should be at least 18 months from delivery date)			
		Sample must be submitted upon participation for accurate evaluation			
		50 Packs (Of 20 Compoules			
6.	Clear Matrix Strip (6mm)	Transparent strips used for composite restorations.			
		Manufactured from mylar polyester film			
		Provides greater strength than acetate			
		Resistant to solvents			
		Does not stretch, curl or break			
		Each strip with a thickness of 0.05mm			
		Sample must be submitted upon participation for accurate evaluation			
		60 Boxes			

**QTN REF: (28) PPN/QTN/2025 (POCD)
 SUPPLY AND DELIVERY OF RESTORATIVE MATERIALS AND INSTRUMENTS FOR COMPOSITE RESTORATION
 PROCEDURES FOR DEPARTMENT OF DENTAL SERVICES, MINISTRY OF HEALTH**

NO	ITEM DESCRIPTION	USER REQUIREMENT	YES	NO	PLEASE SPECIFY THE SPECIFICATION / SCOPE OFFERED FAILURE TO FILL THE OFFER DETAIL SHALL BE IDENTIFIED AS NON-COMPLIANCE TO THE OFFER
7.	Clear Matrix Strip (8mm)	Transparent strips used for composite restorations.			
		Manufactured from mylar polyester film			
		Provides greater strength than acetate			
		Resistant to solvents			
		Does not stretch, curl or break			
		Each strip with a thickness of 0.05mm.			
		Sample must be submitted upon participation for accurate evaluation			
		60 Boxes			
8.	Clear Matrix Strip (10 Mm)	Transparent strips used for composite restorations.			
		Manufactured from mylar polyester film			
		Provides greater strength than acetate			
		Resistant to solvents			
		Does not stretch, curl or break			
		Each strip with a thickness of 0.05mm.			
		Sample must be submitted upon participation for accurate evaluation			
		60 Boxes			