

LAMPIRAN - A1

RUJ: BME/A2K/(1566)/2025/HTD (B4)

1/12/2025

ONE OFF PREVENTIVE MAINTENANCE FOR INTERACOUSTICS AUDIOLOGY EQUIPMENTS FOR AUDIOLOGY UNIT, HOSPITAL RIPAS

SCOPE OF WORK

EQUIPMENT DETAILS (12 UNITS)

BRAND: INTERACOUSTIC

(MODEL, BME NUMBER, SERIAL NUMBER, LOCATION)

1. Affinity 2.0, BME9944, SN:0912358, Audiology unit (RIPASH)
2. Affinity 2.0, BME25574, SN:0991320, Audiology unit (RIPASH)
3. Affinity 2.0, BME25573, SN:0991321, Audiology unit (SSBH)
4. Affinity compact, BME27752, SN:1009014, Child Development Centre (CDC Kiarong Building 3)
5. Equinox PC based Audiometer, BME18621, SN: 0922418, Audiology unit (RIPASH)
6. Equinox PC based Audiometer, BME18620, SN: 0920071, Audiology unit (RIPASH)
7. Titan, BME25994, SN:0911559, Audiology unit (RIPASH)
8. Titan, BME27203, SN:1007324, Child Development Centre (CDC Kiarong Building 3)
9. Titan, BME29938, SN: 1026177, Audiology unit (RIPASH)

10. Eclipse, BME27202, SN:1006383, Audiology unit (SSBH)
11. Eclipse, BME27201, SN:1007765, Child Development Centre (CDC Kiarong Building 3)

SCOPE OF WORK:

Must include but not limited to the following:

- 1) ON-SITE CALIBRATION DONE BY LOCAL CERTIFIED ENGINEER OR OVERSEA ENGINEER WITH FACTORY APPROVED CALIBRATOR
- 2) AIR CONDUCTION
- 3) CALIBRATION OF PURE TONE AIR CONDUCTION
- 4) CALIBRATION OF NARROW BAND MASKING NOISE
- 5) BONE CONDUCTOR CALIBRATION
- 6) CHECKPOINTS:-
 - 6.1) CONTROL DEVICES E.G KNOBS AND PUSH BUTTONS
 - 6.2) INDICATING DEVICES E.G LAMPS AND DISPLAYS CORDS
 - 6.3) CORDS
 - 6.4) PLUGS
 - 6.5) HEADBANDS
 - 6.6) EARPHONES CUSHION
 - 6.7) HOSES AND TUBES FOR IMPEDANCE PROBES
- 7) PROVIDE ALL CERTIFICATION OF CALIBRATION FOR USERS AND BME

LAMPIRAN - A2

RUJ: BME/A2K/(1566)/2025/HTD (B4)

ONE OFF PREVENTIVE MAINTENANCE FOR INTERACOUSTICS AUDIOLOGY EQUIPMENTS FOR AUDIOLOGY UNIT, HOSPITAL RIPAS

SPAREPART

PART NUMBER

QUANTITY