

LAMPIRAN - A1

RUJ: BME/A2K/(1619)/2026/HTD (A17)

25/2/2026

ONE-OFF PPM FOR TWENTY FOUR (24) UNITS OF BLOOD COLLECTION MONITORS IN DEPARTMENT OF LAB SERVICES, MINISTRY OF HEALTH.

SCOPE OF WORK

1. PREVENTIVE MAINTENANCE SHALL FOLLOW EQUIPMENT MANUAL AND BRAND/MODEL SPECIFIC RECOMMENDATIONS AND MAY INCLUDE BUT NOT LIMITED TO:
 - INSPECT THE EXTERIOR OF THE INSTRUMENT
 - INSPECT AND CLEAN THE TRAY
 - ENSURE POWER INLET IS WORKING
 - CHECK MAIN CABLES FOR DAMAGE
 - CHECK THE INSTRUMENT FOR ANY CRACKS IN THE CASEWORK
 - CHECK THE DISPLAY AND KEYBOARD OPERATION
 - PERFORM WEIGHT CALIBRATION USING CALIBRATED STANDARD WEIGHT
2. POST-MAINTENANCE REQUIREMENTS
 - PERFORM FUNCTIONAL CHECKS
 - APPLY PREVENTIVE MAINTENANCE STICKER/TAG AFTER SUCCESSFUL SERVICE
 - DISPOSE OF USED SPAREPARTS/PM KITS UNLESS INSTRUCTED OTHERWISE
3. WORK MUST BE CARRIED OUT BY AUTHORIZED ENGINEERS TRAINED IN MAINTAINING THE EQUIPMENT.
4. VENDORS SHALL SUBMIT A MAINTENANCE CHECKLIST AND SUPPORTING DOCUMENTS SPECIFIC TO THE EQUIPMENT AS PART OF THE TENDER.
5. PREVENTIVE MAINTENANCE MAY BE PERFORMED UP TO ONE MONTH BEFORE THE NEXT SCHEDULED DATE BUT NOT EARLIER.

1	BM	BDC, Bleeding Room	TERUMO	D601	908475	7161	
2	BM	BDC, Bleeding Room	TERUMO	D601	908484	7165	
3	BM	BDC, Bleeding Room	Terumo Penpol	D603	908473	7159	
4	BM	BDC, Bleeding Room	TERUMO	D601	906420	7142	
5	BM	BDC, Bleeding Room	Terumo Penpol	D605	812348	7044	
6	BM	BDC, Bleeding Room	TERUMO	D601	908474	7160	
7	BM	BDC, Mobile	TERUMO	TRAC II	2014 10001	17399	
8	BM	BDC, Bleeding Room	TERUMO	TRAC II	2014 10004	17401	
9	BM	BDC, Bleeding Room	TERUMO	TRAC II	2014 10005	17402	
10	BM	BDC, Bleeding Room	TERUMO	TRAC II	2014 10002	17398	
11	BM	BDC, Mobile	TERUMO	TRAC II	2014 10003	17400	
12	BM	BDC, Bleeding Room	TERUMO	TRAC II	2014 05002	17397	
13	BM	BDC, Bleeding Room	TERUMO	TRAC II	2015 12037	19467	
14	BM	BDC, Bleeding Room	TERUMO	TRAC II	2015 12034	19464	
15	BM	BDC, Mobile	TERUMO	TRAC II	2015 12038	19466	
16	BM	BDC, Mobile	TERUMO	TRAC II	2015 12004	19459	
17	BM	BDC, Bleeding Room	TERUMO	TRAC II	2015 12003	19458	
18	BM	BDC, Mobile	TERUMO	TRAC II	2015 12031	19461	
19	TEMB		TERUMO	D601	908481	7163	
20	TEMB		CENTRON	CM735A	733398	8286	
21	TEMB		CENTRON	CM735A	733399	8285	
22	TEMB		CENTRON	ACCUMIX CM760	76013001	24966	
23	TEMB		CENTRON	ACCUMIX CM760	76013002	24968	
24	TEMB		TERUMO	T-RAC II	2015 12033	19463	

LAMPIRAN - A2

RUJ: BME/A2K/(1619)/2026/HTD (A17)

ONE-OFF PPM FOR TWENTY FOUR (24) UNITS OF BLOOD COLLECTION MONITORS IN DEPARTMENT OF LAB SERVICES, MINISTRY OF HEALTH.

	SPAREPART	PART NUMBER	QUANTITY
1	ONE OFF PPM		24

RUJ: BME/A2K/(1619)/2026/HTD(A17)

ONE-OFF PPM FOR TWENTY FOUR (24) UNITS OF BLOOD COLLECTION MONITORS IN DEPARTMENT OF LAB SERVICES, MINISTRY OF HEALTH

No	District	Location	Brand	Model	Serial No.	No.	Lab Reference No.
1	BRUNEI MUARA	BDC, Bleeding Room	D601		908475	7161	PPM/BDC/BME/2009/002
2	BRUNEI MUARA	BDC, Bleeding Room	TERUMO	D601	908484	7165	PPM/BDC/BME/2009/005
3	BRUNEI MUARA	BDC, Bleeding Room	Terumo Penpol	D603	908473	7159	PPM/BDC/BME/2009/004
4	BRUNEI MUARA	BDC, Bleeding Room	TERUMO	D601	906420	7142	PPM/BDC/BME/2009/007
5	BRUNEI MUARA	BDC, Bleeding Room	Terumo Penpol	D605	812348	7044	PPM/BDC/BME/2009/006
6	BRUNEI MUARA	BDC, Bleeding Room	TERUMO	D601	908474	7160	PPM/BDC/BME/2009/009
7	BRUNEI MUARA	BDC, Mobile	TERUMO	TRAC II	2014 10001	17399	PPM/BDC/BME/2014/047
8	BRUNEI MUARA	BDC, Bleeding Room	TERUMO	TRAC II	2014 10004	17401	PPM/BDC/BME/2014/048
9	BRUNEI MUARA	BDC, Bleeding Room	TERUMO	TRAC II	2014 10005	17402	PPM/BDC/BME/2014/049
10	BRUNEI MUARA	BDC, Bleeding Room	TERUMO	TRAC II	2014 10002	17398	PPM/BDC/BME/2014/050
11	BRUNEI MUARA	BDC, Mobile	TERUMO	TRAC II	2014 10003	17400	PPM/BDC/BME/2014/051
12	BRUNEI MUARA	BDC, Bleeding Room	TERUMO	TRAC II	2014 05002	17397	PPM/BDC/BME/2014/052
13	BRUNEI MUARA	BDC, Bleeding Room	TERUMO	TRAC II	2015 12037	19467	PPM/BDC/BME/2017/067
14	BRUNEI MUARA	BDC, Bleeding Room	TERUMO	TRAC II	2015 12034	19464	PPM/BDC/BME/2017/070
15	BRUNEI MUARA	BDC, Mobile	TERUMO	TRAC II	2015 12038	19466	PPM/BDC/BME/2017/071
16	BRUNEI MUARA	BDC, Mobile	TERUMO	TRAC II	2015 12004	19459	PPM/BDC/BME/2017/072
17	BRUNEI MUARA	BDC, Bleeding Room	TERUMO	TRAC II	2015 12003	19458	PPM/BDC/BME/2017/073
18	BRUNEI MUARA	BDC, Mobile	TERUMO	TRAC II	2015 12031	19461	PPM/BDC/BME/2017/074
19	TEMBURONG		TERUMO	D601	908481	7163	PHM/2011/02
20	TEMBURONG		CENTRON	CM735A	733398	8286	PHM/2010/02
21	TEMBURONG		CENTRON	CM735A	733399	8285	PHM/2012/01
22	TEMBURONG		CENTRON	ACCUMIX CM760	76013001	24966	PPM/PHM/BME/2020/57
23	TEMBURONG		CENTRON	ACCUMIX CM760	76013002	24968	PPM/PHM/BME/2020/56
24	TEMBURONG		TERUMO	T-RAC II	2015 12033	19463	PPM/BB339/17