

**TENDER REFERENCE NO.: KK/76/2026/HTD**

**MINISTRY OF HEALTH  
NEGARA BRUNEI DARUSSALAM**

**TO PPM ONE HUNDRED AND THIRTY NINE (139) UNITS OF  
FRIDGE AND FREEZER IN DEPARTMENT LABORATORY  
SERVICES, MINISTRY OF HEALTH**

**TENDER FEES : \$10.00**

**RECEIPT NO. :**

**CLOSING DATE : 05th May 2026**

**TIME : 2.00 PM**

**FOA :**

**THE CHAIRMAN  
MINI TENDER BOARD, TENDER BOX  
GROUND FLOOR, MINISTRY OF HEALTH  
COMMONWEALTH DRIVE  
BANDAR SERI BEGAWAN BB3910  
NEGARA BRUNEI DARUSSALAM**

**(CLUSTERING)**

**SECTION 2**

**SPECIFICATIONS AND REQUIREMENTS**

**TENDER REFERENCE NO: KK/76/2026/HTD**

**INVITATION TO TENDER**

**ONE-OFF PPM FOR ONE HUNDRED THIRTY-NINE (139) UNITS OF FRIDGES, FREEZERS AND FRIDGE-FREEZER HYBRIDS IN DEPARTMENT OF LAB SERVICES, MINISTRY OF HEALTH**

---

<b>DELIVERY PERIOD</b>	<b>90 DAYS</b>
------------------------	----------------

1. PREVENTIVE MAINTENANCE SHALL FOLLOW EQUIPMENT MANUAL AND BRAND/MODEL SPECIFIC RECOMMENDATIONS AND MAY INCLUDE BUT NOT LIMITED TO:
  - General service
  - Clean and check thermistor
  - Recalibration temperature setting
  - Replace filter and top up gas
  - Inspect cooling coil, compressor and door sealant
  - Clean base plate to ensure to ensure stability
  - Check fan motor and compressor (replace if needed)
  - Check Alarming function (replace battery if needed)
  - Inspect control monitoring display unit
  - Check grease fan motor
2. POST-MAINTENANCE REQUIREMENTS
  - Functional and performance testing
  - Provide a PPM sticker for machines stating company name, time and date of ppm
3. WORK MUST BE CARRIED OUT BY AUTHORIZED ENGINEERS TRAINED IN MAINTAINING THE EQUIPMENT
4. VENDORS SHALL SUBMIT A MAINTENANCE CHECKLIST AND SUPPORTING DOCUMENTS SPECIFIC TO THE EQUIPMENT AS PART OF THE TENDER.
5. PREVENTIVE MAINTENANCE MAY BE PERFORMED UP TO ONE MONTH BEFORE THE NEXT SCHEDULED DATE BUT NOT EARLIER.

No.	District	Location	Brand	Model	Serial No.	BME No.	Lab Reference No.
1	BRUNEI MUARA	BDC	FIOCCHETTI	C65	45452	BME9 815	PPM/BDC/BME/2013/010
2	BRUNEI MUARA	BDC	FIOCCHETTI	C65	45453	BME9 816	PPM/BDC/BME/2013/011
3	BRUNEI MUARA	BDC, CPL	SANYO	MB705 gr	30202309	BME4 999	PPM/BDC/BME/2009/011
4	BRUNEI MUARA	BDC, CPL	HELMER	IPF125	2090033	BME2 5995	PPM/BDC/BME/2022/002
5	BRUNEI MUARA	BDC, CPL	ARCTIKO	ULUF P820 GG	17211107294	BME2 9646	PPM/BDC/BME/2022/003
6	BRUNEI MUARA	BDC, CPL	ESCO	LEXICON II ULT LOW	2012-74729	BME9 695	PPM/BDC/BME/2013/002
7	BRUNEI MUARA	BDC, CPL	SANYO	MDF-U73V VIP SERIES	S/N 50201832	BME5 468	PPM/BDC/BME/2009/010
8	BRUNEI MUARA	BDC, CPL	PHCBI	MDF-DU7OZVX-PE	24120439	BME3 2003	PPM/BDC/BME/2026/001
9	BRUNEI MUARA	BDC, NBTRL	HAIER	HYC-940	BE03L8EIT00QEF5N0004	BME1 8555	PPM/BDC-BTL/15/014
10	BRUNEI MUARA	BDC, NBTRL	HELMER	iBR256-GX	2088053	BME2 6036	PPM/BDC-NBTL/FRI/21/002
11	BRUNEI MUARA	BDC, NBTRL	THERMO SCIENTIFIC	PLR1006	1.01E+13	BME9 904	PPM/BDC-BTL/13/015
12	BRUNEI MUARA	BDC, NBTRL	HAIER	HXC-608	BE06PJT1T00B2E2H0001	BME9 980	PPM/BDC-BTL/14/010
13	BRUNEI MUARA	BDC, NBTRL	SANYO	MDF-U5411	506207181	BME5 644	PPM/BDC-BTL/05/017
14	BRUNEI MUARA	BDC, NBTRL	ESCO	LEXICON UUS-668A-1	2012-63171	BME8 461	PPM/BDC-BTL/12/016

No.	District	Location	Brand	Model	Serial No.	BME No.	Lab Reference No.
15	BRUNEI MUARA	BDC, NBTRL	HELMER	IHB125	1005000	BME8 811	PPM/BDC-BTL/12/011
16	BRUNEI MUARA	BDC, NBTRL	ANTECH	MPR-525	KYC0598GG12507043	BME3 1822	PPM/BDC- NBTL/FRI/25/003
17	BRUNEI MUARA	BDC, NBTRL	HELMER	iBR125-GX	2083438	BME N/A	PPM/BDC- NBTL/FRI/21/001
18	BRUNEI MUARA	BDC, QCL	KIRSCH	BL-V-720	700 8800130	BME6 181	PPM/BDC/BME/2007/001
19	BRUNEI MUARA	BDC, QCL	AUCMA	SC-385E	21157-40 SC-385E 3138-5110-0000-78AV-0044	BME3 1723	PPM/BDC/BME/2022/004
20	BRUNEI MUARA	BDSL	HAIER	HYC-940	BE03L6EIT00B2D4L0006	BME1 4114	PPM/BB165/13
21	BRUNEI MUARA	BDSL	HELMER	IHPF125	091510-46	BME7 684	PPM/BB178/14
22	BRUNEI MUARA	BDSL	HAIER	HYC-940	BE03L100NOOBU99KXJMG	BME1 1234	PPM/BB147/10
23	BRUNEI MUARA	CSRA	BIOBASE	BPR-5V100(G)	123889	BME 27999	PPM/CSRA006/21
24	BRUNEI MUARA	HISTOLOGY LABORATORY	AUCMA	Aucma	BD-99	BME4 722	PPM/HIS/BME/2009/FRIO 02
25	BRUNEI MUARA	HISTOLOGY LABORATORY	KIRSCH	Kirsch	700 88 00115	BME5 645	PPM/HIS/BME/2009/FRI03
26	BRUNEI MUARA	NVRL, MDRUV AUTOCLAVE ROOM	MATRIX	SC-258L	M14262064354000138	BME N/A	PPM/VIRO/2021/031/FF
27	BRUNEI MUARA	MDRUV EXTRACTION ROOM	AUCMA	DW-25W525	5152 5373 9009 2A3W 0035	BME2 5496	PPM/VIRO/2020/023/FF
28	BRUNEI MUARA	NVRL, MDRUV MASTERMIX ROOM	HAIER	DW-30L818	BE0HQ 1E000 0QGM1 J0004	BME N/A	PPM/VIRO/2021/032/FF

No.	District	Location	Brand	Model	Serial No.	BME No.	Lab Reference No.
29	BRUNEI MUARA	NVRL, MDRUV MASTERMIX ROOM	BIOBASE	BPR-541000	YC100025220088	BME2 9670	PPM/VIRO/2022/034/FF
30	BRUNEI MUARA	NVRL, MDRUV OFFICE	ALFA	SC-A186	0073020AIF190917191002126	BME N/A	PPM/VIRO/2020/024/FF
31	BRUNEI MUARA	NCCRL	HAIER	DW-40L508	BE05T7E1T00B2D810018	BME9 886	n/a
32	BRUNEI MUARA	NCCRL	HELMER	IFL 125	SN1004998	BME8 882	PM/CC082/12
33	BRUNEI MUARA	NCCRL	ESCO	LEXICON ULTRA LOW TEMPERATURE	2013-80368	BME9 908	PPM/CC092/13
34	BRUNEI MUARA	NCCRL	HAIER	DW-40L508	BE05T2EAVO00JK9J0014	BME2 4958	PPM/CHE/2020/100
35	BRUNEI MUARA	NCCRL	HAIER	HYCD-282	BEOF18E1TD0B2B880002	BME8 367	PPM/CC081/12
36	BRUNEI MUARA	NCCRL	HAIER	HYC-940	11238	BME1 1238	PPM/CC0075/10
37	BRUNEI MUARA	NCCRL	THERMO SCIENTIFIC	PLR1006	10062013000079	BME9 905	PPM/CC/2019/98
38	BRUNEI MUARA	NCCRL	HAIER	HYC-940	BE03L100N0082 ACBES51	BME7 867	PPM/CC077/11
39	BRUNEI MUARA	NCCRL	ASECO	CHAMPION	30849	BME N/A	PPM/CC033/09
40	BRUNEI MUARA	NCMRL, AMR Room	HAIER	DW-86L288	BE06S9E1T00B2C4S0009	BME9 859	N/A
41	BRUNEI MUARA	NCMRL, Main Laboratory	ASECO	VC200	30858	BME N/A	PPM/MIC/OTH/2009/9
42	BRUNEI MUARA	NCMRL, Main Laboratory	HAIER	HYC-940	BE03L6E1T00B2C1V0004	BME8 823	PPM/MIC/BME/2012/30

No.	District	Location	Brand	Model	Serial No.	BME No.	Lab Reference No.
43	BRUNEI MUARA	NCMRL, Main Laboratory	HAIER	HYC-940	BE03L6E1T00B2C1V0006	BME8 824	PPM/MIC/BME/2012/31
44	BRUNEI MUARA	NCMRL, Serology	HAIER	HYC-610	BE06P 1EOUO OQJL4 Q0002	BME2 4888	PPM/MIC/BME/2020/42
45	BRUNEI MUARA	NCMRL, Serology	MIYAKO	MIS-400S	CK150500042	BME N/A	PPM/MIC/OTH/2016/7
46	BRUNEI MUARA	NCMRL, Serology	HAIER	HYC-1031GD	BE153 4E000 0QJQ3 F0022	BME3 0673	PPM/MIC/BME/2024/54
47	BRUNEI MUARA	NCMRL, Serology	HAIER	DW-40L92	BE02M7EONOOQEECG0005	BME1 8400	PPM/MIC/BME/2015/32
48	BRUNEI MUARA	NHCRLM Andrology Laboratory	THERMO SCIENTIFIC	PLR386	3860091240512	BME9 812	PPM/AND/BME/2018/FRI0 2
49	BRUNEI MUARA	NHCRLM, CYTOLOGY LABORATORY	Thermo Shandon	Thermo Shandon	1006201300084	BME9 903	PPM/CYT/BME/2019/FRI0 5
50	BRUNEI MUARA	NHCRLM, CYTOLOGY LABORATORY	HAIER	HYC-390	BE06M GE1T0 0QEK2 H0003	BME N/A	PPM/CYT/BME/2019/FRI0 8
51	BRUNEI MUARA	NHCRLM, HISTOLOGY	HAIER	HYC-940	N/A	BME1 1236	PPM/HIS/BME/2009/FRI01
52	BRUNEI MUARA	NHCRLM, HISTOLOGY	MIYAKO	MNS450	79000003	BME N/A	PPM/HIS/BME/2009/FRI04
53	BRUNEI MUARA	NHRL	HAIER	Haier HYC-940	BE03L100N00B2AC4QKBH	BME7 869	PPM/HAE-075/1
54	BRUNEI MUARA	NHRL	HAIER	Haier HYC-940	BE03L6E1T00B2C1V0001	BME8 822	PPM/HAE/E0001
55	BRUNEI MUARA	NHRL	HAIER	Haier HYC-940	BE03L8E1T00QEF5N0005	BME1 8556	PPM/HAE/E0038
56	BRUNEI MUARA	NHRL	HAIER	Haier Biomedical	BE153 4E000 OQJQ3 F0021	BME 24960	PPM/HAE/BME/2020/04

No.	District	Location	Brand	Model	Serial No.	BME No.	Lab Reference No.
57	BRUNEI MUARA	NHRL	HAIER	Haier	200242608	BME3 0672	PPM/HAE/BME/2024/01
58	BRUNEI MUARA	NIRL	HAIER	HYC-940	130251578	BME1 1235	PPM/IMM/BME/2009/03
59	BRUNEI MUARA	NIRL	HELMER	HELMER	1004997	BME8 881	PPM/IMM/BME/2005/02
60	BRUNEI MUARA	NIRL	PHCBI	MDF-DU702VH	24060526	BME3 1420	PPM/IMM/BME/2025/01
61	BRUNEI MUARA	NIRL	HAIER	HYCD-282	BE0F18IT00B2BCA001	BME8 368	PPM/IMM/BME/2012/01
62	BRUNEI MUARA	NIRL	THERMO SCIENTIFIC	THERMO SCIENTIFIC CONE000	10062012000004	BME9 727	PPM/IMM/BME/2003/04
63	BRUNEI MUARA	NMRL, SUMBILING	HELMER	ILF120	1003169	BME8 810	PPM/NMRL/BME/003
64	BRUNEI MUARA	NMRL, SUMBILING	HAIER	HYC-610	BE07C1E0N00B2C1G0007	BME8 827	PPM/NMRL/BME/004
65	BRUNEI MUARA	NMRL, SUMBILING	HAIER	HYC-610	BE07C7E1T00B2CCA0003	BME9 849	PPM/NMRL/BME/022
66	BRUNEI MUARA	NMRL, SUMBILING	AUCMA	SC-385E	21157-40SC-385E 3138 5110 0000 78AV 0042	BME N/A	PPM/VIRO/2020/016/FF
67	BRUNEI MUARA	NMRL, SUMBILING	HAIER	HYC-282	BE0F18E1T00QGL3P0003	BME2 4965	PPM/NMRL/BME/029
68	BRUNEI MUARA	NMRL, SUMBILING	LG	GR-V242SL	006INUB00078	BME N/A	PPM/NMRL/BME/025
69	BRUNEI MUARA	NMRL, SUMBILING	AUCMA	BD-220NA	25976-10 BD-220NA 1520 6110 0001 59A7 1248	BME N/A	PPM/NMRL/BME/074
70	BRUNEI MUARA	NVRL EQUIPMENT ROOM	HAIER	HYC-940	BE03L7E1T00QEE9H0004	BME1 7214	PPM/VIRO/2014/015/FF

No.	District	Location	Brand	Model	Serial No.	BME No.	Lab Reference No.
71	BRUNEI MUARA	NVRL EQUIPMENT ROOM, SUMBILING	HAIER	HYC-940	BE03L6E1T00B2C1V0002	BME8 825	PPM/VIRO/2009/011/FF
72	BRUNEI MUARA	NVRL PANTRY ROOM	ALFA		ARM6919080065	BME N/A	PPM/VIRO/2020/029/FF
73	BRUNEI MUARA	NVRL PCR EXTRACTION ROOM	HELMER	ILF125	1004999	BME8 883	PPM/VIRO/2009/012/FF
74	BRUNEI MUARA	NVRL PCR ROOM 1	HAIER	HYC-610	OQJL6CO06	BME2 4889	PPM/VIRO/2020/030/FF
75	BRUNEI MUARA	NVRL ROUTINE LAB ROOM	LG	EXPRESS COOL	5051NKH00613	BME N/A	PPM/VIRO/2009/004/FF
76	BRUNEI MUARA	NVRL SPECIAL TEST ROOM	HAIER	HYC-940	BE03L7E1T00QEE9H0002	BME1 7213	PPM/VIRO/2014/014/FF
77	BRUNEI MUARA	NVRL SPECIAL TEST ROOM	HAIER	HYCD-282	BE0F1 8E1T0 0QGL3 P0006	BME2 4964	PPM/VIRO/2020/019/FF
78	BRUNEI MUARA	NVRL TISSUE CULTURE ROOM	HAIER	HYCD-282	BE0F19E1T00QEE9T0001	BME1 7225	PPM/VIRO/2014/013/FF
79	BRUNEI MUARA	NVRL TISSUE CULTURE ROOM	AUCMA	YC-330	5533 OR73 9009 2A37 0618	BME2 5446	PPM/VIRO/2020/022/FF
80	BRUNEI MUARA	NVRL TISSUE CULTURE ROOM	ARCTIKO	LAF 700-40C / DAI 0290-4	1114737	BME2 9648	PPM/VIRO/2022/033/FF
81	BRUNEI MUARA	NVRL, REAGENT PREPARATION ROOM	SANYO	MDF-436	50811315	BME5 854	PPM/VIRO/2009/001/FF
82	BRUNEI MUARA	NVRL, MCRUV RECEPTION ROOM	AUCMA	BD-220NA	25976-10 BD-220NA 1520 6110 0001 59A7 1248	BME N/A	PPM/VIRO/2020/017/FF
83	BRUNEI MUARA	NVRL, MDRUV EXTRACTION ROOM	HAIER	DW-86L828J	BE0F0 0EAW0 0QGL2 MOOO1	BME2 4961	PPM/VIRO/2020/020/FF
84	BRUNEI MUARA	NVRL, MDRUV MASTERMIX ROOM	AUCMA	DW-25W525	5152 5373 9009 2A3W 0034	BME2 5497	PPM/VIRO/2020/018/FF

No.	District	Location	Brand	Model	Serial No.	BME No.	Lab Reference No.
85	BRUNEI MUARA	NVRL,SPECIAL TEST ROOM , SUMBILING	HITACHI	R-570AG5	2150800078	BME5 851	PPM/VIRO/2009/003/FF
86	BRUNEI MUARA	POCT	HAIER	HYC -282	BE0F1 8E1TO OQGLE P0001	BME2 4963	PPM/POC/RF/01
87	BRUNEI MUARA	QAPSU	HAIER	HYC-68	BE069 0E000 0QEQ6 70026	BME3 1393	PPM/CAL/2025/BME/01
88	BRUNEI MUARA	TB LAB	AUCMA	BD-400NA	1537 5010 0001 59AE 0490	BME N/A	PPM/VIRO/2020/025/FF
89	TUTONG	Assay 2 Suite, MGS	BIOBASE	BPR-5V310	YCO31025212252	BME2 7865	PPM/MGS/2022/025
90	TUTONG	Assay 2 Suite, MGS	BIOBASE	BDF-25V268	DF25V26825210160A	BME2 7862	PPM/MGS/2022/022
91	TUTONG	Assay 3 Suite, MGS	BIOBASE	BDF-25V268	DF25V26825210161A	BME2 7863	PPM/MGS/2022/023
92	TUTONG	Assay 1 Suite, MGS	BIOBASE	BDF-25V268	DF25V26825210162A	BME2 7864	PPM/MGS/2022/024
93	TUTONG	Assay 3 Suite, MGS	BIOBASE	BPR-5V310	YCO31025211763	BME2 7866	PPM/MGS/2022/026
94	TUTONG	Assay 1 Suite, MGS	BIOBASE	BPR-5V310	YCO31025211764	BME2 7867	PPM/MGS/2022/027
95	TUTONG	PMMH Micro Lab	HELMER	HLF120	2073454	BME2 4969	PPM/PMM/BME/20/P001
96	TUTONG	PMMH Blood bank	HAIER	DW-86L628E	BE06Q OEAYO OQGM3 50030	BME2 7545	-
97	TUTONG	PMMH Micro Lab	HAIER	HYC-940	BE03L6E1T00B2C1V0008	BME8 826	PPM/PMM/BME/12/P003
98	TUTONG	Main Lab, MGS	ARCTIKO	ULUF P820 GG	1106973	BME2 9645	PPM/MGS/BME/2022/1

No.	District	Location	Brand	Model	Serial No.	BME No.	Lab Reference No.
99	TUTONG	Main Lab, MGS	ARCTIKO	LAF 700	1114736	BME2 9649	PPM/MGS/BME/2022/2
100	TUTONG	Main Lab, MGS	HAIER	HYC-610	BE07C 2EAQ0 0QJM3 J0012	BME2 7167	PPM/MGS/BME/2022/4
101	TUTONG	Main Lab, MGS	FIOCCHETTI	VISION 2T 700 ECT-F TOUCH (VST0702)	82706	BME2 9708	PPM/MGS/BME/2022/3
102	TUTONG	BSL3, CMDLID	ARCTIKO	LF100	1117109	BME2 9644	PPM/NMRL/BME/063
103	TUTONG	BSL2, CMDLID	ARCTIKO	ULUF P820 GG	35211112867	BME2 9647	PPM/NMRL/BME/067
104	TUTONG	SPECIMEN RECEIVING, CMDLID	FIOCCHETTI	VST 0702	82704	BME2 9706	PPM/NMRL/BME/069
105	TUTONG	BSL2, CMDLID	FIOCCHETTI	VST 0702	82705	BME2 9707	PPM/NMRL/BME/070
106	TUTONG	BSL2, CMDLID	HAIER	HYC-610	BE07C 2EAQ0 0QJM3 J003	BME2 7164	PPM/NMRL/BME/049
107	TUTONG	BSL3, CMDLID	HAIER	HYC-610	BE07C 2EAQ0 0QJM3 H0026	BME2 7165	PPM/NMRL/BME/050
108	TUTONG	BSL3, CMDLID	HAIER	HYC-610	BE07C 2EAQ0 0QJM3 J004	BME2 7166	PPM/NMRL/BME/051
109	TUTONG	BSL3, CMDLID	HAIER	HYC-610	BE07C 2EAQ0 0QJM3 J0005	BME2 7168	PPM/NMRL/BME/053
110	TUTONG	PMM Blood Bank	REVCO	Reb-23040	N14G-3234	BME2 0748	PPM/PMM/BME/97/P017
111	TUTONG	PMM Blood Bank	KIRSCH	Super 720	7000110547	BME2 0747	PPM/PMM/09/P009
112	TUTONG	PMM Clinical Chem Lab	KIRSCH	Labex-332	3300201090	BME2 0746	PPM/PMM/97/P018

No.	District	Location	Brand	Model	Serial No.	BME No.	Lab Reference No.
113	TUTONG	PMM Clinical Chem Lab	KIRSCH	Labex	3301000373	BME5 784	PPM/PMM/25/P002
114	TUTONG	PMM Clinical Chem Lab	HAIER	HYC-940	BE03L 100N0 OBU99 US6SU	BME1 1237	PPM/PMM/09/P010
115	TUTONG	PMM Micro Lab	HAIER	HYC-940	BE03L7E1T00B2D4L00001	BME2 0745	PPM/PMM/13/P011
116	TUTONG	PMM Micro Lab	HAIER	HYC-940	BE03L7E1 T00QEE9H0001	BME1 7259	PPM/PMM/14/P014
117	TUTONG	PMM Micro Lab	HAIER	HYC-941	BE03L7E1TOOQEE9H0005	BME1 7260	PPM/PMM/14/P015
118	BELAIT	SSBH, MICROBIOLOGY	HAIER	HYC-940	BE032100N00B2ACH2FOE	BME2 0592	PPM/SSB/2011/09
119	BELAIT	SSBH, MICROBIOLOGY	ANTECH	MPR-525	KYC0598G12507028	BME3 1807	PPM/SSB/2025/36
120	BELAIT	SSBH, CLINICAL CHEMISTRY	HAIER	HYC-910	BE03L100N00BU99UTT6M	BME1 1233	PPM/SSB/2009/03
121	BELAIT	SSBH, CLINICAL CHEMISTRY	HAIER	HYC-940	BE03L NE0N0 0QJL3 L0001	BME2 3040	PPM/SSB/2020/83
122	BELAIT	SSBH, CLINICAL CHEMISTRY	MATRIX	M2SC-850X	WM22067-1	BME N/A	PPM/SSB/2023/106
123	BELAIT	SSBH, CLINICAL CHEMISTRY	ANTECH	MPR-525	KYC0598G12507026	BME3 1805	PPM/SSB/2025/37
124	BELAIT	SSBH, HAEMATOLOGY	HAIER	HYC-910	BE07C2E0N00B2D7D0001	BME9 888	PPM/SSB/2013/07
125	BELAIT	SSBH, HAEMATOLOGY	HAIER	HYC-910	BE07C2E0N00B2D7D0008	BME9 889	PPM/SSB/2013/08
126	BELAIT	SSBH, BLOOD TRANSFUSION	SANYO	1404GR	2006-16192	BME1 1253	PPM/SSB/2009/06

No.	District	Location	Brand	Model	Serial No.	BME No.	Lab Reference No.
127	BELAIT	SSBH, BLOOD TRANSFUSION	HAIER	DW-40L508	BE05T7E1TOOB20810013	BME9 887	PPM/SSB/2013/14
128	BELAIT	SSBH, BLOOD TRANSFUSION	ANTECH	MPR-525	KYC0598G12507027	BME3 1806	PPM/SSB/2025/38
129	BELAIT	SSBH BLOOD TRANSFUSION	ANTECH	MPR-525	KYCO598612411009	BME3 1007	PPM/SSB/2025/04
130	BELAIT	SSBH, SAMPLE PROCESSING (SEPARATION)	THERMO SCIENTIFIC	PL6500	3860089240512	BME9 729	PPM/SSB/2013/02
131	BELAIT	SSBH, SAMPLE PROCESSING (SEPARATION)	AUCMA	PL6500	36130-10 BD-702 (1C60 6110 0000 3B3N 0352)	BME2 6970	
132	BELAIT	SSBH, BLOOD BANK	HELMER	HBR125-GX	2161194	BME3 1419	PPM/SSB/2025/20
133	BELAIT	SSBH, BLOOD BANK	HELMER	HBR125-GX	2004831	BME9 852	PPM/SSB/2013/10
134	BELAIT	SSBH, BLOOD BANK	EXLOLEXICA N ULT			BME2 0334	PPM/SSB/2013/15
135	TEMBURONG	PIHMH	HAIER	HYCD-282	BE0F18E1T00QEE850003	BME1 8401	PPM/PHM/2015/40
136	TEMBURONG	PIHMH	KIRSCH	SUPER-720	7000110576	BME2 0704	PPM/PHM/2009/20
137	TEMBURONG	PIHMH	HELMER	IB125	2002649	BME9 562	PPM/PHM/2013/34
138	TEMBURONG	PIHMH	HELMER	HBR125-GX	2161193	BME3 1418	PPM/PHM/2025/03
139	TEMBURONG	PIHMH	HAIER	HYCD-282	BE0F18E1T00QGL2B0001	BME 23181	PPM/PHM/BME/2020/55

**SECTION 3**  
**TENDER FORM**

To:

**TENDER REFERENCE NO: KK/76/2026/HTD**

**INVITATION TO TENDER**  
**ONE-OFF PPM FOR ONE HUNDRED THIRTY-NINE (139) UNITS OF FRIDGES, FREEZERS AND**  
**FRIDGE-FREEZER HYBRIDS IN DEPARTMENT OF LAB SERVICES, MINISTRY OF HEALTH**

<b>DELIVERY PERIOD</b>	
<b>SCOPE OF WORK</b>	<b>PRICE (B\$)</b>
<ol style="list-style-type: none"><li>1. PREVENTIVE MAINTENANCE SHALL FOLLOW EQUIPMENT MANUAL AND BRAND/MODEL SPECIFIC RECOMMENDATIONS AND MAY INCLUDE BUT NOT LIMITED TO:<ul style="list-style-type: none"><li>▪ General service</li><li>▪ Clean and check thermistor</li><li>▪ Recalibration temperature setting</li><li>▪ Replace filter and top up gas</li><li>▪ Inspect cooling coil, compressor and door sealant</li><li>▪ Clean base plate to ensure to ensure stability</li><li>▪ Check fan motor and compressor (replace if needed)</li><li>▪ Check Alarming function (replace battery if needed)</li><li>▪ Inspect control monitoring display unit</li><li>▪ Check grease fan motor</li></ul></li><li>2. POST-MAINTENANCE REQUIREMENTS<ul style="list-style-type: none"><li>▪ Functional and performance testing</li><li>▪ Provide a PPM sticker for machines stating company name, time and date of ppm</li></ul></li><li>3. WORK MUST BE CARRIED OUT BY AUTHORIZED ENGINEERS TRAINED IN MAINTAINING THE EQUIPMENT</li><li>4. VENDORS SHALL SUBMIT A MAINTENANCE CHECKLIST AND SUPPORTING DOCUMENTS SPECIFIC TO THE EQUIPMENT AS PART OF THE TENDER.</li><li>5. PREVENTIVE MAINTENANCE MAY BE PERFORMED UP TO ONE MONTH BEFORE THE NEXT SCHEDULED DATE BUT NOT EARLIER.</li></ol>	
<b>TOTAL PRICE (B\$)</b>	

NO.	SPAREPART	PART NUMBER	QUANTITY
1	ONE OFF PPM		139