

LAMPIRAN - A1

RUJ: BME/A2K/(1682)/2026/HTD (B9)

11/5/2026

ONE-OFF PREVENTIVE MAINTENANCE FOR ONE (1) UNIT OF DENTAL CHAIR COMPLETE SYSTEM, MODEL: COMPACT i5, S/N: UPCV574500 AT ISTANA NURUL IZZAH CLINIC

SCOPE OF WORK

Vendor must be an authorised agent to provide maintenance services.

1. Perform Air Compressor inspection including checking regulated air pressure and draining excess condensation water.

2. Perform Metasys suction inspection including cleaning air filter and cleaning any clogged waste.

3. Perform full preventive maintenance replacements for the dental chair which include:

- Oring
- the return air filter
- white filter
- sealing rubber
- the valve body with the viton O ring
- foot control case assembly
- the control membrane
- the pedal knob
- the rubber foot
- water filter membrane
- Durr disposable filter
- air filter membrane
- burket spray water valve
- bowl filter
- cover cap
- vohran tube for diameters 10 mm and 16.5mm.

4. Replace the green and clean WK solution bottle with a new one.

5. Perform inspection for mounted Pro X of the dental chair

6. Perform general checking of the dental chair after preventive maintenance is done. This includes checking;

- water pressure
- air pressure and suction
- water and air inlet valves
- syringes
- low/high speed turbine handpiece
- scalers
- headrest motor
- display monitor
- and overall condition of all the handpieces and their hoses.

7. Perform general cleaning of the dental chair by using disinfectant wipes to clean the patient seat and the general surface area.

LAMPIRAN - A2

RUJ: BME/A2K/(1682)/2026/HTD (B9)

ONE-OFF PREVENTIVE MAINTENANCE FOR ONE (1) UNIT OF DENTAL CHAIR COMPLETE SYSTEM, MODEL: COMPACT i5, S/N: UPCV574500 AT ISTANA NURUL IZZAH CLINIC

	SPAREPART	PART NUMBER	QUANTITY
1	Maintenance Kit Compact i5 WEK	30016974	1
2	Maintenance Kit Planmeca Compact i5	30016585	1
3	Green and clean Wk Solution (4-Pack)	10035933	1
4	Perform Metasys WEK service kit replacement		