

BIL	Quotation Reference	Description	Advertisement Date	Closing Date (Not Later Than 09.00AM)	Quotation Fee	Requesting Department	Focal Person
2	45/QTN/JPK/2024	<p>SUPPLY AND DELIVERY OF TWO HUNDRED (200) PIECES OF ADULT PEAK FLOW METER FOR ALL HEALTH CENTRES, DEPARTMENT OF HEALTH SERVICES</p> <p><u>PLACE OF SUBMISSION:</u> QUOTATION BOX (GROUND FLOOR) MINISTRY OF HEALTH COMMONWEALTH DRIVE BANDAR SERI BEGAWAN, BB 3910 NEGARA BRUNEI DARUSSALAM</p>	04/03/2025	22/03/2025	\$5.00	JABATAN PERKHIDMATAN KESIHATAN, KEMENTERIAN KESIHATAN.	HAJAH NORHALIZAWATY BINTI HJ ABD RAZAK PEGAWAI JURURAWAT TEL: 2382652

TAWARAN SEBUTHARGA
 JABATAN PERKHIDMATAN KESIHATAN

Rujukan Tawaran: **45/QTN/JPK/2024**

Nama Tawaran:

SUPPLY AND DELIVERY OF TWO HUNDRED (200) PIECES OF ADULT PEAK FLOW METER FOR ALL HEALTH CENTRES, DEPARTMENT OF HEALTH SERVICES

SPECIFICATIONS	
STANDARD FEATURES	<p>ADULT PEAK FLOW METER - 200 PIECES A peak flow meter helps monitor respiratory conditions such as asthma, by measuring maximum airflow out of lungs. This device used during consultation with doctor or asthma nurse educators.</p> <p>SPECIFICATIONS</p> <ul style="list-style-type: none"> - To provide accurate peak expiratory flow (PEF) measurements as part of asthma management plan - Suitable for clinic use - Suitable for use in children and adult - Available with EU scale (tested to the European standard) and/or ATS scale (tested to the American Thoracic Society standard) - Measurement Range: approximately 60 L/min to minimum of 800 L/min (EU scale or ATS scale) - Accuracy: approximately $\pm 10\%$ or 10 L/min (EU Scale) or $\pm 10\%$ or 20 L/min (ATS Scale) - Each device is individually calibrated for absolute accuracy - Easy-to-read scale (text on a different colour background) displays in scale increments of 10 L/min until 700 L/min - Portable and lightweight - Durable - All Latex free materials - Built-in one-way valve preventing accidental inhalation through the device. - Can be cleaned for multi-patient use - Complete with a detachable, reusable and sterilisable mouthpiece - Suitable for multiple patient use with availability of consumables to minimize cross infection between patients such as Universal Sterilisable Mouthpieces, Adult Disposable Expiratory One Way Valve Mouthpieces, Two Way Mouthpieces

• Authorized by:

Dr Sia Ai Tee

TAWARAN SEBUTHARGA
JABATAN PERKHIDMATAN KESIHATAN

Rujukan Tawaran: **45/QTN/JPK/2024**

Nama Tawaran:

SUPPLY AND DELIVERY OF TWO HUNDRED (200) PIECES OF ADULT PEAK FLOW METER FOR ALL HEALTH CENTRES, DEPARTMENT OF HEALTH SERVICES

SPECIFICATIONS	
STANDARD FEATURES	TWO HUNDRED (200) PIECES OF ADULT PEAK FLOW METER FOR ALL HEALTH CENTRES, DEPARTMENT OF HEALTH SERVICES

Please state:

- Brand/ Model
- Accessories
- Country of origin
- International Safety Standard
- Year Manufactured
- Shelf Life
- Price per unit
- Total price
- Delivery period
- Price Validity
- Sample & Brochure
- Warranty

• Authorized by:

Dr Sia Ai Tee

**TAWARAN SEBUTHARGA
JABATAN PERKHIDMATAN KESIHATAN**

Rujukan Tawaran: **45/QTN/JPK/2024**

Nama Tawaran:

**SUPPLY AND DELIVERY OF TWO HUNDRED (200) PIECES OF ADULT PEAK FLOW METER FOR
ALL HEALTH CENTRES, DEPARTMENT OF HEALTH SERVICES**

SPECIFICATIONS		YES	NO	IF 'NO' PLEASE STATE ALTERNATIVE OFFER
STANDARD FEATURES	To provide accurate peak expiratory flow (PEF) measurements as part of asthma management plan			
	Suitable for clinic use			
	Suitable for use in children and adult			
	Available with EU scale (tested to the European standard) and/or ATS scale (tested to the American Thoracic Society standard)			
	Measurement Range: approximately 60 L/min to minimum of 800 L/min (EU scale or ATS scale)			
	Accuracy: approximately $\pm 10\%$ or 10 L/min (EU Scale) or $\pm 10\%$ or 20 L/min (ATS Scale)			
	Each device is individually calibrated for absolute accuracy			
	Easy-to-read scale (text on a different colour background) displays in scale increments of 10 L/min until 700 L/min			
	Portable and lightweight			
	Durable			
All Latex free materials				

STANDARD FEATURES	Built-in one-way valve preventing accidental inhalation through the device			
	Can be cleaned for multi-patient use			
	Complete with a detachable, reusable and sterilisable mouthpiece			
	Suitable for multiple patient use with availability of consumables to minimize cross infection between patients such as Universal Sterilisable Mouthpieces, Adult Disposable Expiratory One Way Valve Mouthpieces, Two Way Mouthpieces			

PROCUREMENT REQUIREMENT	
Brand/ Model	
Country Of origin	
International Safety Standard	
Year Manufactured	
Shelf Life (minimun 1.5 Years)	
Price per unit	
Total price	
Delivery period	
Price Validity	
Sample & Brochure	
Warranty	