

Rujukan Kami: (87) MOH/HQ/P/IKLAN-SH/2025

LAMPIRAN 3

BIL	Quotation Reference	Description	Advertisement Date	Closing Date (Not Later Than 09.00AM)	Quotation Fee	Requesting Department
3	BME/A2K/(1449)/2025/HTD (A20)	REPAIR ONE (1) UNIT OF INDUCTIVELY COUPLED PLASMA MASS SPECTROMETRY SYSTEM (ICPMS) BRAND: AGILENT, MODEL: 7700, S/N: JP15163366, BME18768 AT ENVIRONMENT LABORATORY, DEPARTMENT OF SCIENTIFIC SERVICES.	11/06/2025	05/07/2025	\$5.00	JABATAN TEKNOLOGI PENJAGAAN KESIHATAN, KEMENTERIAN KESIHATAN

## **LAMPIRAN - A1**

**RUJ: BME/A2K/(1449)/2025/HTD (A20)**

**3/6/2025**

**REPAIR ONE (1) UNIT OF INDUCTIVELY COUPLED PLASMA MASS SPECTROMETRY SYSTEM (ICPMS) BRAND: AGILENT, MODEL: 7700, S/N: JP15163366, BME18768 AT ENVIRONMENT LABORATORY, DEPARTMENT OF SCIENTIFIC SERVICES.**

### **SCOPE OF WORK**

Repair Agilent 7700 ICPMS at Environment Laboratory, Department of Scientific Service, inclusive of but not limited to:

- Scope of work to follow equipment manual/recommendation specific for this system which include the following where if necessary:
  - (i) installation of items listed in spareparts list (Lampiran A2)
  - (ii) alignment
  - (iii) calibration
  - (iv) optimisation
  - (vi) any other related maintenance works required after installation of items listed in spareparts list
- Perform functional test, verification and performance check post-maintenance.
- Disposal of any used spare parts, unless otherwise instructed.
- Work on the system is to be carried out by engineer(s) authorised in providing maintenance to the equipment.

## LAMPIRAN - A2

RUJ: BME/A2K/(1449)/2025/HTD (A20)

REPAIR ONE (1) UNIT OF INDUCTIVELY COUPLED PLASMA MASS SPECTROMETRY SYSTEM (ICPMS) BRAND: AGILENT, MODEL: 7700, S/N: JP15163366, BME18768 AT ENVIRONMENT LABORATORY, DEPARTMENT OF SCIENTIFIC SERVICES.

SPAREPART		PART NUMBER	QUANTITY
1	Argon flow controller, 7700x/s (Pb-Free)	G3280-67135	1 unit