

Rujukan Kami: (96) MOH/HQ/P/IKLAN-SH/2025

LAMPIRAN 3

BIL	Quotation Reference	Description	Advertisement Date	Closing Date (Not Later Than 09.00AM)	Quotation Fee	Requesting Department
3	BME/A2K/(1473)/2025/HTD (B8)	TO REPAIR ONE (1) UNIT OF HEAT SEALER (STERILE POUCH), BRAND: GETINGE, MODEL: 6003000195 AT CSSD RIPAS HOSPITAL.	08/07/2025	26/07/2025	\$5.00	JABATAN TEKNOLOGI PENJAGAAN KESIHATAN, KEMENTERIAN KESIHATAN.

LAMPIRAN - A1

RUJ: BME/A2K/(1473)/2025/HTD (B8)

26/6/2025

TO REPAIR ONE (1) UNIT OF HEAT SEALER (STERILE POUCH), BRAND: GETINGE, MODEL: 6003000195 AT CSSD RIPAS HOSPITAL.

SCOPE OF WORK

To repair one (1) unit of heat sealer (sterile pouch), Brand: Getinge, Model: 6003000195 at CSSD Ripas Hospital. Detail are as follow:

1. BME26018, Brand: Getinge, Model: 6003000195, SN: PPR20210056, at CSSD, W/order 2025-6832

- Scope of work to follow equipment manual/recommendation specific for this system which include the following where if necessary:

- (i) installation of items listed in spareparts list (Lampiran A2)
- (ii) alignment
- (iii) calibration
- (iv) optimisation
- (vi) any other related maintenance works required after installation of items listed in spareparts list (Lampiran A2)

- Perform functional test, verification and performance check post-maintenance.

- Disposal of any used spare parts, unless otherwise instructed.

- Work on the system is to be carried out by engineer(s) authorised in providing maintenance to the equipment.

LAMPIRAN - A2

RUJ: BME/A2K/(1473)/2025/HTD (B8)

**TO REPAIR ONE (1) UNIT OF HEAT SEALER (STERILE POUCH), BRAND: GETINGE, MODEL: 6003000195 AT
CSSD RIPAS HOSPITAL.**

SPAREPART		PART NUMBER	QUANTITY
1	PRINTER MODULE PROSEAL	68000103	1

WORK ORDER

WORKORDER NO: 2025-6832		BME NO: BME26018	
USER WORKORDER NO: HTD/20202		JOB TYPE: <input checked="" type="checkbox"/> CM <input type="checkbox"/> PM <input type="checkbox"/> IIC <input type="checkbox"/> OB <input type="checkbox"/> TR	
USERNAME: Abdullahled		DATE REPORTED: 14/06/2025 TIME:	
TEL: 6490		EQUIPMENT: HEAT SEALER (STERILE POUCH)	
LOCATION: CSSD		MANUFACTURER: GETINGE	
BME: RIPAS HOSPITAL		MODEL: 6003000195	
PRIORITY: <input checked="" type="checkbox"/> NORMAL <input type="checkbox"/> CRITICAL		SERIAL NO: PPR20210056	
ASSIGNED BY: TO:		DATE ASSIGNED: TIME:	
REPORTED PROBLEM: CANNOT PRINT			
TECHNICAL REPORT	SUPPLIER: MC QUIPP TECHNISSELL		Date Started: Time Started:
OBSERVED PROBLEM:		ACTION TAKEN: Contacted local supplier for further action and check 26/6/2025 - Rcvd qtn \$3,250.00 Proceed Iklan	
ACTION SUMMARY: <input checked="" type="checkbox"/> Waiting for inspection report / quotation <input type="checkbox"/> Process obsolete <input type="checkbox"/> Inspection Report (Signed by user) <input type="checkbox"/> Equipment taken by supplier <input checked="" type="checkbox"/> Quotation received [Date: 26 June 2025] [Date:] <input checked="" type="checkbox"/> Process B2K / A2K / IKLAN			
REF: ()BME/A2K/HTD/() ()BME/B2K/HTD/() BME/A2K()/HTD()		OB. REF: ()BME/OBS/HTD/()	
JOB STATUS: <input type="radio"/> OPEN <input type="radio"/> PENDING <input checked="" type="radio"/> IN PROGRESS <input type="radio"/> COMPLETE <input type="radio"/> OBSOLETE		Date: Time:	
Engineer's Name and Signature :		User's Name, Signature & Stamp :	

Date :