

LAMPIRAN - A1

RUJ: BME/A2K/(1526)/2025/HTD (B8)

1/10/2025

TO REPAIR FIVE (5) UNITS OF CRITICAL BED, BRAND: HILLROM, MODEL: TC-155 & P1900 AT CRITICAL CARE UNIT RIPAS HOSPITAL.

SCOPE OF WORK

To repair five (5) units of Critical bed, Brand: hillrom, Model: TC-155 & P1900 at Critical care unit Ripas Hospital. Detail are as follow:

1. BME6480, Brand: Hillrom, Model: P1900, SN: 1234AM3528, at ICU 2, W/order 2025-8033 (Item 1)
2. BME18790, Brand: Hillrom, Model: P1900, SN: Q030AM0024, at ICU 2, W/order 2025-7683 (Item 2)
3. BME5481, Brand: Hillrom, Model: TC-155, SN: G066AM5742, at ICU 2, W/order 2025-7153 (Item 1,3,4,5)
4. BME3485, Brand: Hillrom, Model: 1900A, SN: B018AM1302, at ICU 1, W/order 2025-7385 (Item 1,6,7)
5. BME5482, Brand: Hillrom, Model: P1900, SN: G066AM5744, at ICU 2, W/order 2025-6149 & 2024-1895 (Item 1,6,7,8)

- Scope of work to follow equipment manual/recommendation specific for this system which include the following where if necessary:

- (i) installation of items listed in spareparts list (Lampiran A2)
- (ii) alignment
- (iii) calibration
- (iv) optimisation
- (vi) any other related maintenance works required after installation of items listed in spareparts list (Lampiran A2)

- Perform functional test, verification and performance check post-maintenance.

- Disposal of any used spare parts, unless otherwise instructed.

- Work on the system is to be carried out by engineer(s) authorised in providing maintenance to the equipment.

LAMPIRAN - A2

RUJ: BME/A2K/(1526)/2025/HTD (B8)

TO REPAIR FIVE (5) UNITS OF CRITICAL BED, BRAND: HILLROM, MODEL: TC-155 & P1900 AT CRITICAL CARE UNIT RIPAS HOSPITAL.

	SPAREPART	PART NUMBER	QUANTITY
1	VLV SOLENOID 2W/2P DELTROL (P1900)	6543035	8
2	POWER CONTROL PCB ASSY, 230V/50HZ TC-155	130067	1
3	RH CAREGVR BD & GASKET-TC	SA1732	1
4	LEFT SIDE RAIL INT. CABLE	SA1735	1
5	RIGHT SIDE RAIL INT. CABLE	SA1734	1
6	RESERVOIR W/SEALS DELTROL	127940	2
7	HYDRAULIC FLUID-CHART	36199S	4
8	EMERGENCY TRENDELENBURG VALVE	6543038	1