

LAMPIRAN - A1

RUJ: BME/A2K/(1665)/2026/HTD (A12)

21/4/2026

TO PERFORM PLANNED PREVENTIVE MAINTENANCE SERVICING FOR MINDRAY BENEHEART D3 DEFIBRILLATORS, BENEHEART D1 AEDs AND BENEHEART C1A AEDs IN MINISTRY OF HEALTH

SCOPE OF WORK

Scope of work on the Defibrillator and Automatic External Defibrillator (AED) to follow Mindray manual/recommendation specific for model D3 defibrillators, Beneheart C1A Fully Automatic & Beneheart D1, which include the following where necessary:

- (i) inspection
- (ii) cleaning
- (iii) functional test
- (iv) performance test
- (v) battery calibration
- (vi) any other related preventive maintenance works required.

- Functional check post-maintenance.

- Apply preventive maintenance sticker/tag after successful service.

- Work on the system is to be carried out by engineer(s) authorised in providing maintenance to Mindray Defibrillator.

- Vendor(s) to submit maintenance checklist/details/scope specific for Mindray Defibrillator along with other necessary documents to join this tender.

By-

LAMPIRAN - A2

RUJ: BME/A2K/(1665)/2026/HTD (A12)

TO PERFORM PLANNED PREVENTIVE MAINTENANCE SERVICING FOR MINDRAY BENEHEART D3 DEFIBRILLATORS, BENEHEART D1 AEDs AND BENEHEART C1A AEDs IN MINISTRY OF HEALTH

	SPAREPART	PART NUMBER	QUANTITY
1	PPM FOR BENEHEART D3		13
2	PPM FOR D1 AED		22
3	PPM FOR C1A AED		31