

LAMPIRAN - A1

RUJ: BME/A2K/(1692)/2026/HTD (B4)

2/6/2026

TO PERFORM PLANNED PREVENTIVE MAINTENANCE SERVICING FOR MINDRAY N15 PATIENT MONITOR, N12 PATIENT MONITOR , N1 PATIENT MONITOR, TM80 TELEMETRY MONITOR AND BENEVISION CENTRAL MONITOR STATION AND WORKSTATION AT CCU,RIPASH

SCOPE OF WORK

Scope of work on the Patient Monitor/Telemetry Monitor/Central Monitor to follow Mindray manual/recommendation specific for model N15/N12/N1 Patient Monitor, TM80 Telemetry Monitor and Benevision Central Monitor, which include the following where necessary:

- (i) inspection
 - (ii) cleaning
 - (iii) functional test
 - (iv) performance test
 - (v) battery inspection
 - (vi) 24/7 Service call for 1 year Period
- any other related preventive maintenance works required.

- Functional check post-maintenance.
- Apply preventive maintenance sticker/tag after successful service.
- Work on the system is to be carried out by engineer(s) authorised in providing maintenance to Mindray Patient Monitor/Telemetry Monitor/Central Monitor
- Vendor(s) to submit maintenance checklist/details/scope specific for Mindray Patient Monitor/Telemetry Monitor/Central Monitor along with other necessary documents to join this tender.

LAMPIRAN - A2

RUJ: BME/A2K/(1692)/2026/HTD (B4)

TO PERFORM PLANNED PREVENTIVE MAINTENANCE SERVICING FOR MINDRAY N15 PATIENT MONITOR, N12 PATIENT MONITOR , N1 PATIENT MONITOR, TM80 TELEMETRY MONITOR AND BENEVISION CENTRAL MONITOR STATION AND WORKSTATION AT CCU,RIPASH

	SPAREPART	PART NUMBER	QUANTITY
1	PM for N15		20
2	PM For N12		24
3	PM for N1		44
4	PM for TM80		20
5	PM FOR HP(MINDRAY) CMS		2
6	PM FOR HP(MINDRAY) REMOTE CENTRAL VIEWSTATION		4
7	COMPREHENSIVE MAINTENANCE		1

LAMPIRAN – A3

RUJ: BME/A2K/(1692)/2026/HTD (B4)

TO PERFORM PLANNED PREVENTIVE MAINTENANCE SERVICING FOR MINDRAY N15 PATIENT MONITOR , N12 PATIENT MONITOR , N1 PATIENT MONITOR, TM80 TELEMETRY MONITOR AND BENEVISION CENTRAL MONITOR STATION AND WORKSTATION AT CCU, RIPASH

NO	ITEM DESCRIPTION	PART NUMBER	QUANTITY
1	Brand: Mindray® Equipment: Patient Monitor Model: N15 Manufacturer: Shenzhen, China Location: RIPASH CCU Serial Nos: 1. F5-34055163/BME30430 2. F5-34055164/BME30431 3. F5-34055165/BME30433 4. F5-34055166/BME30433 5. F5-34055167/BME30434 6. F5-34055168/BME30435 7. F5-34055169/BME30436 8. F5-34055170/BME30437 9. F5-34055171/BME30438 10. F5-34055172/BME30439		
2	Brand: Mindray® Equipment: Patient Monitor Model: N12 Manufacturer: Shenzhen, China Location: RIPASH CCU Serial Nos: 1. F8-34072264/BME30440 2. F8-34072265/BME30441 3. F8-34072266/BME30442 4. F8-34072267/BME30443 5. F8-34072268/BME30444 6. F8-34072268/BME30445 7. F8-34072270/BME30446 8. F8-34072271/BME30447 9. F8-34072272/BME30448 10. F8-34072273/BME30449 11. F8-34072274/BME30450 12. F8-34072275/BME30451		

3	Brand: Mindray® Equipment: Patient Monitor Model: N1 Manufacturer: Shenzhen, China Location: RIPASH CCU Serial Nos: 1. AA3-34072703/BME30470 2. AA3-34072704/BME30471 3. AA3-34072705/BME30472 4. AA3-34072706/BME30473 5. AA3-34072707/BME30474 6. AA3-34072708/BME30475 7. AA3-34072709/BME30476 8. AA3-34072710/BME30477 9. AA3-34072711/BME30478 10. AA3-34072712/BME30479 11. AA3-34072713/BME30480 12. AA3-34072714/BME30481 13. AA3-34072715/BME30482 14. AA3-34072716/BME30483 15. AA3-34072717/BME30484 16. AA3-34072718/BME30485 17. AA3-34072719/BME30486 18. AA3-34072720/BME30487 19. AA3-34072721/BME30488 20. AA3-34072722/BME30489 21. AA3-34072723/BME30490 22. AA3-34072724/BME30491		
4	Brand: Mindray® Equipment: Telemetry Monitor Model: TM80 Manufacturer: Shenzhen, China Location: RIPASH CCU Serial Nos: 1. HJ-35018074/BME30460 2. HJ-35018068/BME30461 3. HJ-35018071/BME30462 4. HJ-35018077/BME30463 5. HJ-35018069/BME30464 6. HJ-35018070/BME30465 7. HJ-35018076/BME30466 8. HJ-35018072/BME30467 9. HJ-35018073/BME30468 10. HJ-35018075/BME30469		

5	Brand: Mindray® Equipment: Central Monitor Station and Workstation Model: Benevision CHP-mindray) Manufacturer: Shenzhen, China Location: RIPASH CCU Serial Nos: 1. AA3-34072728/BME30452 2. AA3-34072729/BME30457 3. AA3-34072730/BME30458		
6	Comprehensive Maintenance of CMS, Workstation and Datamed PC. 1 Lot Replace any defective part of Central Monitoring System, Workstation and Datamed (12 Lead ECG Converter) *Comprehensive Maintenance offer is only valid for 1 year PPM period proposal. *		